



Blue™ is no longer available for purchase.
Contact us to order Steam Chex™ today!
customerservice@proppermfg.com
Tel. (718) 392-6650

Reorder No.: Quantity:
26810800 48 Rolls • ¾" x 60 YDS
26810900 36 Rolls • 1" x 60 YDS

Blue Autoclave Indicator Tape and Steam Chex are pressure-sensitive process indicator tapes designed to seal packs that are exposed to steam sterilization.

Blue Autoclave Indicator Tape and Steam Chex are made of saturated crepe paper and coated with solvent-free adhesive.

Blue Autoclave Indicator Tape and Steam Chex are not made with lead based ink or natural latex rubber.

Blue Autoclave Indicator Tape features a **WHITE** to dark brown/black color change
Steam Chex Indicator Tape features a **GREEN** to dark brown/black color change

