

Propper Laryngoscopes are designed for Longevity, Ease-of-Use and Infection Control!

The Laryngoscope System is

- **Versatile** : 3 handle sizes available with two light sources and 11 blade varieties.
- **Easy To Sterilize** : Handles feature removable inner compartments allowing hassle free cleaning and sterilization. Blades are designed with integrated fiber optic bundles and made of matte polished select quality stainless steel to eliminate crevices and prevent fluid entrapment.
- **Standardized** : Blades and Handles are Green System compliant.

Stubby Handle				“C” Handle				“AA” Handle			
											
Hook Head	Outer Shell	Battery Compartment	Light Switch	Hook Head	Outer Shell	Battery Compartment	Bottom Cap	Hook Head	Outer Shell	Battery Compartment	Bottom Cap
Reorder No.				Reorder No.				Reorder No.			
19958100	19957800			19958200	19915000			19958300	19915200		
Stubby LED	Stubby Xenon			C LED	C Xenon			AA LED	AA Xenon		

Laryngoscope Handles

- Standard Battery Sizes
- Xenon or Superlume LED bulbs
- Ball and lock system and removable inner compartments prevent fluid entrapment allowing hassle free disinfection and sterilization
- Compatible with Steam, EO and Vaporized Hydrogen Peroxide sterilization
- Fits any blade compatible with “Green” ISO 7376 standard

Superlume® LED

- Brightest output (400% brighter than Halogen)
- High color temperature for increased visibility
- LED light source lasts 1,000x longer than a Halogen bulb
- LED’s long lifetime will make frequent bulb replacements a thing of the past

Propper Laryngoscope Blades

Propper Manufacturing blades feature an integrated, high-performance fiber optic bundle. The large diameter fiber optics, combined with state-of-the-art Superlume LED handles, provides an unmatched visualization for laryngoscope users.

Propper Laryngoscope blades are made of select quality 304 Stainless Steel with flawless pinhole free surface finish. The epoxy sealed integrated fiber optic channel is heat resistant to 392°F. This streamlined system minimizes fluid entrapment and enables sterilization in Steam, EO and Vaporized Hydrogen Peroxide.

Propper blades feature an easy clip system and are Green System compliant.

Propper Laryngoscope Blades' light output will not fall below 2,000 lux for five years when used with Propper LED handles.

Macintosh Blade



Miller Blade



Size	Description	Reorder No.
1	Macintosh Blade	19916100
2	Macintosh Blade	19914400
3	Macintosh Blade	19914500
4	Macintosh Blade	19914600
5	Macintosh Blade	19916200
0	Miller Blade	19915600
1	Miller Blade	19915700
1.5	Miller Blade	19915500
2	Miller Blade	19915800
3	Miller Blade	19915900
4	Miller Blade	19916000

Disclaimer: The laryngoscope blades are matte polished and handles are plated to improve the surface quality. Deep etching, mechanical impacts, excessive force used in manual cleaning, especially from dried biological materials, deviations from the recommended point-of-use care, disinfection and sterilization procedures may result in damage to the instrument surface.

