

Propper Manufacturing Company has been supporting you with best-in-class sterilization assurance products and services for nearly a century.

Now, we proudly introduce

STERI-WRAP® **Unite**™

SMS bonded wrap for use in steam, ethylene oxide and vaporized hydrogen peroxide sterilization cycles.

Steri-Wrap Unite 60 is available in sheets

31502400	24" x 24"	1 Case of 120
31503600	36" x 36"	1 Case of 72
31504800	48" x 48"	1 Case of 48
31505400	54" x 54"	1 Case of 48



718-392-6650 
customerservice@proppermfg.com 
www.proppermfg.com 

STERI-WRAP® Unite™

Propper Manufacturing Company Steri-Wrap Unite is a dual color, bonded SMS sterilization wrap for easy-to-use, simultaneous wrapping.

Benefits and Characteristics

- ∧ Barrier allows sterilant penetration but prevents penetration of microorganisms
- ∧ Double-wrap protection in half the time compared with sequential double wrapping
- ∧ Strong, tear resistant 60 gsm material soft enough to enable easy wrapping
- ∧ Dual color bonded blue/green wrap helps differentiate specific instrument sets
- ∧ Lot code printed on every application
- ∧ FDA approved and CE marked
- ∧ Use in steam, ethylene oxide and hydrogen peroxide sterilization
- ∧ Shelf Life After Sterilization: 365 days

STERI-WRAP® **Unite™** 60



Reorder No.

31502400

120 Applications/Case
24in x 24in | 61cm x 61cm

31503600

72 Applications/Case
36in x 36in | 91cm x 91cm

31504800

48 Applications/Case
48in x 48in | 122cm x 122cm

31505400

48 Applications/Case
54in x 54in | 137cm x 137cm