



**propper**  
Manufacturing Co., Inc.

Contact us to order Verline today!  
[customerservice@proppermfg.com](mailto:customerservice@proppermfg.com)  
Tel. (718) 392-6650

# Introducing **VERLINE™**

A pressure-sensitive process indicator  
tape designed to seal packs that are  
exposed to steam sterilization.

Reorder No.:

Quantity:

**26810500**

**48 Rolls • ¾" x 60 YDS**

**26810600**

**36 Rolls • 1" x 60 YDS**



@proppermfg